



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Surgical Mobile Light

Product Code :
NMOTL10006



Description :

Surgical Mobile Light Doma Dia: 20in, 90,000 Lux:
Halogen Bulb: 24V - 150W
Reflector: Diachronic
Colour Temp.: 4300K
Temp. Increase at Head Area: 2-4C
Finish: Oven Baked

Surgical Mobile Light Doma Dia: 14in, 50,000 Lux:
Halogen Bulb: 24V - 100W
Reflector: Diachronic
Colour Temp.: 4300K
Temp. Increase at Head Area: 2-4C
Finish: Oven Baked Surgical Mobile Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

