

Product Name :
Face Mask With SS Ring**Product Code :**
NMANI10017**Description :**

Face Mask With SS Ring Size 3. Black Rubber. Padded. CE certified

Technical Specification :

Anesthesia mask. Medical mouth piece mask. Rubber. 5 point metal ring for easy and tight fixation with head straps (not included).

Ideal for use with all types of automatic and manual resuscitators or for anaesthesia.

Ergonomically modelled for a perfect fit on the patient's face and optimized grip for the operator.

With inflatable cushion to ensure perfect fit around the mouth.

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,