

Product Name :
Reservoir Bag, PVC with Adapter

Product Code :
NMANI10031



Description :

Reservoir Bag, PVC with Adapter To fit our ambu type bags codes: AN002c, AN002f, AN002g & AN002h

Technical Specification :

Non-Rebreather Mask with Reservoir-Adult Size
Low-resistance check valves prevent rebreathing and allow exhaled gas to escape
Safety vent allows entrainment of room air
Clear, soft vinyl for patient comfort and visual patient assessment

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)