



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :

Ambu Type Bag (Artificial Resuscitator), Adult
(Autoclavable) Black Rubber

Product Code :

NMANI10036

**Description :**

Ambu Type Bag (Artificial Resuscitator), Adult (Autoclavable) Black Rubber Autoclavable. Optional : Reservoir bag - AN403 (Order separately). CE certified

Technical Specification :

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.



Medical Hospital Equipments Manufacturers

Website: SITE