



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**

Ambu Type Bag (Artificial Resuscitator), (Autoclavable)  
Superior, Adult, Black Rubber

**Product Code :**

NMANI10039



**Description :**

Ambu Type Bag (Artificial Resuscitator), (Autoclavable) Superior, Adult, Black Rubber Black rubber. Optional :  
Reservoir bag - AN403 (Order separately). CE certified

**Technical Specification :**

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.



**Medical Hospital Equipments Manufacturers**

**Website:** Medical Hospital Equipments Manufacturers,