



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ambu Type Bag (Artificial Resuscitator), (Autoclavable)
Superior, Adult, Black Rubber

Product Code :
NMANI10039



Description :

Ambu Type Bag (Artificial Resuscitator), (Autoclavable) Superior, Adult, Black Rubber Black rubber. Optional : Reservoir bag - AN403 (Order separately). CE certified

Technical Specification :

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.



Medical Hospital Equipments Manufacturers

Website: SITE