



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :

Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable

Product Code :

NMANI10041

**Description :**

Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable Black rubber. With face mask size. 2. Round shape. CE certified

Technical Specification :

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,