



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Corrugated Tube - Black Rubber

**Product Code :**  
NMANI10043



**Description :**

Corrugated Tube - Black Rubber Size: 42

**Technical Specification :**

Material: Rubber  
Colour: Black  
Excellent flexibility and wearability, high-strength, flame-retardant.  
Commonly used for special machines.



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,