



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Head Harness Rubber

**Product Code :**  
NMANI10045



**Description :**

Head Harness Rubber Several holes permit ventilation and holds both the Child & Adult heads. It consists of four rubber straps for attaching to the metal Mask Ring, which fits over the facemask

**Technical Specification :**

It is designed to hold both Child & Adult heads.

It covers a wide range of applications.

Range includes supplied-air respirators, air-purifying respirators.

Powered air-purifying respirators.



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,