

**Product Name :**  
Head Harness Silicon**Product Code :**  
NMANI10046**Description :**

Head Harness Silicon Color available : Black, White

**Technical Specification :**

Silicone Head Harness features center vents for ventilation and reinforced holes for durability.

Silicone Head Strap is non-sterile and molded for comfort.

Silicone Face Mask Harness consists of four straps which provide good support of face mask and evenly distributes pressure over a large area.

With minimal spacing between the holes on the straps, these designs allow for a wide range of locations to secure the straps to any sized individual.

Silicone Head Harness is not made with natural rubber latex.

