

Product Name :
Anesthesia machine with LED Screen

Product Code :
NMANI10054



Description :

Anesthesia machine with LED Screen With LED screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers . Handle for easy maneuvering Epoxy painted high quality mild steel frame. Equipped with Anesthesia Vaporizer. Agent type: Isoflurane, Enflurane, Halothane. Classical CO2 circle absorber with 1 canister With necessary monitoring and alarm function, to ensure the operational safety

Technical Specification :

Suitable for adults and pediatric patients
Equipped with Large LCD screen, various monitoring waveform and parameters indifferent colors
Pneumatically driven electronically controlled
Integrated CO2 absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast release APL valve
6-tube of flow meters for 3-gas sources meet long operation requirement
Spirometry loops
Equipped with Anesthesia Vaporizer.

