

Product Name :
Anesthesia Machines with CO2, Multi Parameter,
Anesthesia Depth Monitor

Product Code :
NMANI10056



Description :

Anesthesia Machines with CO2, Multi Parameter, Anesthesia Depth Monitor With 5.5 SNT screen. Mounted on 4 castors. 2 front castors with brakes. With 1 safety lock and 2 big drawers Handle for easy maneuvering. Aluminium alloy rail frame. With durable and non-corrosive stainless steel platform. Equipped with high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise. The agent type: Isoflurane, Enflurane. With independent Anesthesia Ventilator, the ventilator frame is moulded and height can be adjusted. The ventilator can display P-t and V-t wave-form at the same time. Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage. Fresh gas outlet is in front of frame for the user to easily connect Jackson-Rees circuit. With upgraded functions, eg, CO2 Monitor, Multi-parameter monitor, Anesthesia depth monitor, etc. Optional : O2, N2O Yoke device. When main gas supply is fails, the backup gas supply starts.

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,