

Product Name :
Anesthesia Machine with 10.4 TFT screen

Product Code :
NMANI10057



Description :

Anesthesia Machine with 10.4 TFT screen With 10.4 TFT color screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers. Handle for easy maneuvering. Epoxy painted high quality mild steel frame. Equipped with Two high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise. The agent type: Isoflurane, Enflurane, Halothane, Sevoflurane. Classical CO₂ circle absorber with 2 canisters for low resistance and minimum gas leakage. Optional: O₂, N₂O Yoke device. When main gas supply is fails, the Backup gas supply starts.

Technical Specification :

Electromagnetic proportional valve control technology and precision flow sensor monitoring technology. Precisely pressure controlled ventilation mode and automatically adjusted flow rate enable it provides rapid gas leakage compensation for all kinds of patients. Automatic compensation maintains accurate tidal volumes, and ensures that what you set is what you get. Comprehensive respiratory mechanical parameters monitor, real-time waveforms display, optimizing clinical decision. Intelligent alarm management. Electronic PEEP function.

