



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Ceiling Surgical Operating Light Small Single

Product Code :
NMOTL10012



Description :

Ceiling Surgical Operating Light Small Single has daylight characters and can be adjusted as per the operation requirements. Lights arms are made of Stainless Steel Fitting. The range of the color temp - 4,200 K. Detachable handle for simplistic sterilization. Parabolic reflectors & color correcting filters.

Dome Head: 550 mm

Halogen Bulbs: 24V x 70 W x 4 nos. Light intensity at 1 mt: 50,000 Lux

Height Adjustment: 914 mm

Parabolic Reflector: 180 mm

Color Temperature: 4200 K

Control Unit: Transformer Ceiling Surgical Operating Light Small Single, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE