

**Product Name :**  
Rack Steam Sterilizer, Single Rack**Product Code :**  
NMALV10016**Description :**

Rack Steam Sterilizer, Single Rack Steam sterilizer, pressure type, single rack, used to sterilize reusable: syringes and needles. Aluminium vessel with high pressure seal suitable for carrying out sterilisation using superheated steam. Steam sterilizer designed to operate on external heat source (non-electric) at  $+121\frac{1}{2}^{\circ}\text{C}$  over a 20 minute period. Components: Aluminium vessel with heat resistant handles. Aluminium cover with heat resistant handles, rubber sealing gasket, pressure valve & safety plug. 1 aluminium rack with 1 lid, for 36 syringes ( 4 syringes 5ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles. Sterilization capacity: 5.5 litres  
Supplied with: 1 x sterilization indicator sticker (pack of 100), 2 x spare safety plugs, 3 x spare sealing gaskets, 1 x plastic bowl for washing syringes, 1 x carrying bag, 1 x user manual in English.

**Technical Specification :**

Made of Aluminium

For placing of syringes inside autoclave

Holds 36 syringes (4 syringes 5 ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles

With Lid

