

**Product Name :**  
Rack Steam Sterilizer, Double Rack**Product Code :**  
NMALV10017**Description :**

Rack Steam Sterilizer, Double Rack Steam sterilizer, pressure type, double rack, used to sterilize reusable: syringes and needles. Aluminium vessel with high pressure seal suitable for carrying out sterilization using superheated steam. Steam sterilizer designed to operate on external heat source (non-electric) at  $+121\frac{1}{2}^{\circ}\text{C}$  over a 20 minute period. Components: Aluminium vessel with handles. Aluminium cover with handles and rubber sealing gasket. Pressure valve. Safety plug. 2 aluminium racks with 2 lid, for 64 syringes (44 syringes 2 ml, 16 syringes 1 ml, 4 syringes 5 ml) and 80 needles per rack. Sterilization capacity: 9 litres. Supplied with: 1 x sterilization indicator sticker (pack of 100), 2 x spare safety plugs, 3 x spare sealing gaskets, 1 x plastic bowl for washing syringes, 1 x carrying bag, 1 x user manual in English.

**Technical Specification :**

Made of Aluminium

For placing of syringes inside autoclave

Holds 36 syringes (4 syringes 5 ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles

With Lid

