

Product Name :
Accessories for Rack Sterilizer**Product Code :**
NMALV10018**Description :**

Accessories for Rack Sterilizer

Technical Specification :

Three different types of stainless steel containers of different heights (large, medium and small) are available for handling different loads.

Stacking shorter containers can increase load capacity (two medium-sized and three small-sized containers can fit in a given chamber).

A mesh wire basket is suitable for most loads, including glassware and media.

A pail or a mesh wire basket with a solid bottom is used for agar and to contain potential spillage from a torn biohazard bag.

A drum has a hinged lid with a handle. Depending on the load, it serves as a closed (but not watertight) trash container or as a carrying and storage case.

The side and bottom of the drum are perforated for better steam penetration during sterilization.

Sliding panels cover these perforations prior to and after sterilization.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,