

Product Name :
Autoclave 10 Litre Electric**Product Code :**
NMALV10022**Description :**

Autoclave 10 Litre Electric Stainless Steel 304 seamless construction. Auto cut-off heating element and indication light. With time adjustment facility. Dimension : 235mm Dia x 235mm Ht. (with Lid) Note: Height and capacity adjustable as per customers requirement.

Technical Specification :

Autoclaves are sturdily constructed in rustproof, heavy gauge aluminum for use in rugged condition.

Its a unique vessel that generates steam under pressure and temperature between 120°C to 135°C.

On reaching the right temperature the instruments are subject to a saturated steam environment for a pre-specified time which results in safe and guaranteed sterilization.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)