

Product Name :
Autoclave 10 Litre Non-Electric

Product Code :
NMALV10023



Description :

Autoclave 10 Litre Non-Electric Stainless Steel 304 seamless construction. Dimension : 235mm Dia x 240mm Ht. (with Lid) Note: Height and capacity adjustable as per customers requirement

Technical Specification :

Autoclaves are sturdily constructed in rustproof, heavy gauge aluminum for use in rugged condition.

Its a unique vessel that generates steam under pressure and temperature between 120°C to 135°C.

On reaching the right temperature the instruments are subject to a saturated steam environment for a pre-specified time which results in safe and guaranteed sterilization.

NAUGRA®

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,