



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Autoclave 11 Litre Non Electric

Product Code :
NMALV10025



Description :

Autoclave 11 Litre Non Electric Stainless Steel 304 seamless construction. Dimension : 220mm Dia x 300mm Ht
Note: Height and capacity adjustable as per customers requirement

Technical Specification :

Type: Pressure Cooker Type/Laboratory bench autoclave(N-Type). Dimensions: 9" Dia X 11" Height .
Equipped with pressure gauge, Steam Release Valve, safety valve, pressure regulation device, with heat resistance lid handle and body handle, inner stand.
Ideal Use: To sterilize medical, equipment, glassware, surgical, equipment, medical wastes, culture media, autoclave containers, plastic tubes, and pipette tips.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,