

Product Name :
Stainless Steel Wire Mesh Trays And Wire Baskets**Product Code :**
NMALV10027**Description :**

Stainless Steel Wire Mesh Trays And Wire Baskets These sterilizing Trays & baskets made in Stainless Steel 304, allow a systematic organization of instrument for sterilization. The comfortable handling provides safety and control through all fields of applications in sterilization. Storage and transportation are also perfectly organized with these baskets and trays Available in various shapes, sizes & designs to accommodate wide range of Surgical Instruments & Equipment

Technical Specification :

Stainless Steel Wire Mesh Baskets for your specific application requirements.

Available in different mesh openings, construction materials, and widths, these are widely used in kitchens and medical fields as colander and disinfection basket respectively.

As all the wire ends are melted to a spherical shape, there is no risk of cut, scratch or graze while handling.