



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Cidex Tray, Stainless Steel

**Product Code :**  
NMALV10029



**Description :**

Cidex Tray, Stainless Steel Size : 705x204x150mm

**Technical Specification :**

Cidex Tray is manufactured with stainless steel 202 grade with thickness of 0.50mm without having any sharp edges.

These trays have seamless jointless construction.

Cidex trays are used to accomodate Medical Instruments and Endoscopes.

Stainless Steel Cidex Tray is durable and Reusable.



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,