

Product Name :
Electric - S.S. 304 Grade**Product Code :**
NMALV10038**Description :**

Electric - S.S. 304 Grade Instrument Sterilizers, electric, with hinged lid fitted with immersion rod and water outlet tap, seamless stainless steel construction. Fitted with automatic cutoff safety pin device, thermostat & tray lifting arrangement

Technical Specification :

Instrument Sterilizer is made up of seamless stainless steel of High Quality 304 Grade.
Sterilizer has Automatic tray Lifting system.
Handle on side with insulated grip for proper safety.
The Instrument Sterilizer is thermostatically controlled.
The Instrument Sterilizer uses 220/240 V, 50 Cycles, Single Phase A.C.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)