

**Product Name :**  
Hospital Dressing Drums - Stainless Steel**Product Code :**  
NMALV10043**Description :**

Hospital Dressing Drums - Stainless Steel Dressing Drum, domed top and bottom with sliding removable band and drop down-hasps (stainless steel)

**Technical Specification :**

Dressing drums at VMS are available in different sizes  
Products that have been mentioned as Super Deluxe are made up of heavy quality  
High-quality material is used for engineering  
Highly reliable and durable  
Can be used to store sterile bathrobes, dressings, and operational linen etc.