

Product Name :
Hospital Dressing Drums - Stainless Steel**Product Code :**
NMALV10044**Description :**

Hospital Dressing Drums - Stainless Steel Dressing Drum, domed top and bottom with sliding removable band and drop down-hasps (stainless steel).

Technical Specification :

Dressing drums at VMS are available in different sizes
Products that have been mentioned as Super Deluxe are made up of heavy quality
High-quality material is used for engineering
Highly reliable and durable
Can be used to store sterile bathrobes, dressings, and operational linen etc.