

Product Name :
Laboratory Rectangular Heating Plate with Cast Iron Top**Product Code :**
NMALV10049**Description :**

Laboratory Rectangular Heating Plate with Cast Iron Top This is very useful for continuous heating upto 350 \pm 1/2 $^{\circ}$ C. The body is made of specially designed with perforated P.C.R.C. sheet painted with attractive stove enamel.

Technical Specification :

These hot plates are manufactured using high grade raw material and consists of cast iron table top, frame made of P.C.R.C. sheet coated powder paint and energy controller to regulate the temperature from ambient up to 350 degree C.

Moreover, the table top of the unit consists of heating elements of 1.5 kW rating and is round in shape with diameter of 20 cm.