

**Product Name :**  
Laboratory Autoclaves - Semi Automatic Series

**Product Code :**  
NMALV10051



**Description :**

Laboratory Autoclaves - Semi Automatic Series Microprocessor based PID controller with dual display for Temp. & Time Unique pedal free design Reduced loading height No hinge or shaking lids Fully SS 304 construction 30 PSI models also available on request.

**Technical Specification :**

Processor based Controller with time and temperature programmable by a user.  
Operating and Display range up to 122oC.  
Pt100 sensor, enabling precise monitoring.  
Low water level alarm and cut off.  
Sensor open alarm.  
End of cycle buzzer and auto reset.  
Miniature Circuit Breaker for On/Off plus safety.  
Solenoid Valve for automatic purging and exhaust in a fully automatic