



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infant Weighing Scale - Beam Type

Product Code :
HW20



Description :

Infant Weighing Scale - Beam Type Beam type (mechanical) scale for infants, 16 kg x 10 gm. Measuring range up to approx 16 kg. Minimum graduation, 10 g. Readout in kg only. Easy readable in low light working situations. Reading time max 15 seconds. Includes removable tray for infant. With easy zero adjustment possibility. Design allows rough handling. Smooth surface/finishing allows for easy cleaning/disinfection. All vital parts made of rust proof materials. Easy maintenance and repair in low tech settings. Splash proof and shock resistant light-weight body. Materials; body, metal (white baked enamel) or plastic Sturdy design allows heavy duty use.



Medical Hospital Equipments Manufacturers

Website: [SITE](#)