



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Hot Water Bottles

Product Code :
MRG20100006



Description :

Capacity: 1 ltr, 2 ltr, 2.5 ltr, 3 ltr

100% leak proof.

Seamless moulded.

Made from first grade Natural Rubber only.

Fitted with heat resistant Polypropylene screw-stopper with a rubber safety seal for preventing water leakage.

Also available in metal stopper.

Easy and smooth filling.

Very supple and soft feel. Hot Water Bottles, Medical Rubber Goods Exporters, Medical Rubber Goods Manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)