



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Shadowless OT Light with Single Dome Dia 521mm

Product Code :
NMOTL10019



Description :

Ceiling Shadowless OT Light with Single Dome Dia 521mm The Light Head provides the superior performance and is good for all kinds of surgeries. Light arms are made of Stainless Steel Fitting The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Sterilizable handle has been provided on the dome for easy movement Heat absorbing & color correcting filter glass. Technical specifications
Dome head : 521mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 1014mm Color Temperature : 4200 K \pm 300 Temp rise in field : 2° - 3° C Control Unit : Transformer



Medical Hospital Equipments Manufacturers

Website: SITE