

**Product Name :**  
High Vacuum Suction Machine SS

**Product Code :**  
NMSUC10HV2



#### Description :

Complete S.S Body to prevent from rust. Heavy Duty ½ hp motorized pump for fast vacuum built – up. Absolutely Noiseless. Self lubricating system for trouble free service. Finely regulable 0-730mm Hg vacuum, Easy to read off. Over Flow Safety Device: over – flow check valve and removable filter to protect the pump from moisture. Non collapsible tubing.

Max Vacuum: 730 ± 10mm Hg regulable through Vacuum Control Knob.

Suction Capacity: 40 to 45 ltrs. / MIN

Vacuum Gauge\Bourden Type 4” 0-760mm Hg caliberation

Suction Bottle: Glass jars or 2 x 2000cc Poly carbonate jars ( UNBREAK ABLE ) with over flow check valve (float type)

Motor: 1/2 hp motor

Vacuum Pump: Rotary type.

Power Source: AC 220V 50 Hz. Single Phase.

Overall Size: L 380mm x W 390mm x H 825mm (Approx.)

Finish: Powder coated.

Weight: 34 kg ( appox ) High Vacuum Suction Machine SS, High Vacuum Suction Exporters, High Vacuum Suction Manufacturers

---

**Website:**Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**