

Product Name :
Led Shadowless Operation Lamp**Product Code :**
NMOTL10050**Description :**

Led Shadowless Operation Lamp Technical Specifications: Illumination : 160000 LUX Color temperature : 4300±500K Color indication index (Ra) : 95 Power of LED bead : N128x1W Spot diameter d50 : 130mm-350mm Illumination depth : 140cm Illumination-adjustment range : 10%~100% Ultraviolet irradiance with wavelength under 400nm : 0.002W/m² Focusing control : Manual Perfect for endoscopy surgery Power supply : 220V, 50Hz /110V, 60Hz LED average service life : 60000hrs Lowest height of installation : 2800mm Available in ceiling, floor and wall configurations

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,