



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Led Shadowless Operation Lamp

Product Code :
NMOTL10051



Description :

Led Shadowless Operation Lamp Spot center illumination : 120000LUX Color temperature : 4300±500K Color indication index (Ra) : 95 Power of LED bead : N64x1W Spot diameter d50 : 130mm-350mm Illumination depth : 120cm Illumination-adjustment range : 10%~100% Ultraviolet irradiance with wavelength under 400nm : 0.002W/m² Focusing control : manual Power supply : 220V, 50Hz /110V, 60Hz LED average service life : 60000hrs Lowest height of installation : 2800mm Available in ceiling, floor and wall configurations



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,