

Product Name :
Led Surgical Light Mobile - 3 Lamp**Product Code :**
NMOTL10052**Description :**

Led Surgical Light Mobile - 3 Lamp Diameter of lamp head : Main lamp, -600mm Illumination : Main lamp, 70000 lx ~ 160000 lx, Subsidiary Color temperature : 3000k ~ 6000k Color rendering index : 92 ~ 100 Diameter of Spot: D10 Main lamp 120mm, Subsidiary : D50 Main lamp 70mm, Subsidiary Depth of illumination : = 600mm Adjustment of Height : 800mm Brightness control : 9 steps Total radiation illumination : = 500W/M2 Power supply : A.C100 ~ 220V 50 HZ Power consumption : 50VA Power of bulb : 3 WLED Number of LED Bulbs : Main lamp, 21

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com