

Product Name :
High Vacuum Suction LIPO**Product Code :**
NMSUC10HV3**Description :**

Complete S.S Body to prevent from rust. Heavy Duty 1 hp motorized pump for fast vacuum built – up. Absolutely Noiseless. Self lubricating system for trouble free service. Finely regulable 0-750mm Hg vacuum, Easy to read off. Over Flow Safety Device: over – flow check valve and removable filter to protect the pump from moisture. Non collapsible tubing

Max Vacuum: 730 ± 10mm Hg regulable through Vacuum Control Knob.

Suction Capacity: 40 to 45Lt per Min

Vacuum Gauge Bourden Type 4 0-750mm Hg calibration

Suction Bottle: Glass jars or 2 x 2000cc Poly carbonate jars (UNBREAK ABLE) with over flow check valve (float type)

Motor: 1 hp motor

Vacuum Pump: Rotary type (double pump) .

Power Source: AC 220V 50 Hz. Single Phase.

Overall Size: L 380mm x W 390mm x H 825mm (Approx.)

Finish: Powder coated.

Weight: 34 kg (appox) High Vacuum Suction LIPO, High Vacuum Suction Exporters, High Vacuum Suction Manufacturers

Website:Medical Hospital Equipments Manufacturers,

NAUGRA[®]

Medical Hospital Equipments Manufacturers