



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Suction Tube

**Product Code :**  
NMSUC100P2



**Description :**

Suction Tube: Silicon or PVC transparent tubing with soft sealing grips. Retain elasticity even after prolonged use. ID =10mm & wall thickness = 2.5 mm. Excellent flexible & heat resistant. Odor less & repeatedly autoclaveable. Transparency of tube allows visibility of flowing fluids Suction Tube, Parts Suction Units Exporters, Parts Suction Units Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,