



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mobile Shadowless Operation Lamp Single Dome Stand
Model With Castors & With Battery Backup (4 Bulbs)

Product Code :
NMOTL10074



Description :

Mobile Shadowless Operation Lamp Single Dome Stand Model With Castors & With Battery Backup (4 Bulbs)
Illumination Intensity : 40000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height
Adjustment : 400mm Lamphead Diameter : 500mm Power Supply : 220V / 50Hz Power Input : 120 VA With
battery backup



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,