



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
LED Examination Light (Wall Mounted)

Product Code :
NMOTL10082



Description :

LED Examination Light (Wall Mounted) Illumination Intensity : 12000 Lux Color Temperature : 6000-6500 K Spot Diameter : 80 - 250mm Power Supply : 220 V / 50 Hz. AC Extention : 75 cm goose neck Specification of Bulb : LED 1x3W Brightness Control : Available



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,