

Product Name :
LED Examination Light With Stand Floor Model & With
Microcomputer Controlled 8 Step Brightness Adjustment

Product Code :
NMOTL10083



Description :

LED Examination Light With Stand Floor Model & With Microcomputer Controlled 8 Step Brightness Adjustment
Standard Features : Microcomputer control 8-steps brightness adjustment Technical specifications Color
Temperature : 3000-5500 K Illumination : ? 20000 Lux Temperature Rise of Field : ?8(°C) Power Supply : 220 V
/ 50 Hz. AC Power (W) : ?12 Type of bulb : Super bright LED Luminotron 3V/0.5Wx19 Spot Diameter : 80-120
mm Range : 300-500mm from lamp adapter Brightness Control : 8 step digital control

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,