

**Product Name :**  
LED Examination Light with Stand Floor Model &  
Microcomputer Controlled 8 Step Brightness Adjustment  
With Battery Backup

**Product Code :**  
NMOTL10084



**Description :**

LED Examination Light with Stand Floor Model & Microcomputer Controlled 8 Step Brightness Adjustment With Battery Backup Standard Features: Microcomputer control 8-steps brightness adjustment Technical specifications Color Temperature : 3000-5500 K Illumination : ? 20000 Lux Temperature Rise of Field : ?8(°C) Power Supply : 220 V / 50 Hz. AC Power (W) : ?12 Type of bulb : Super bright LED Luminotron 3V/0.5Wx19 Spot Diameter : 80-120 mm Range : 300-500mm from lamp adapter Brightness Control : 8 step digital control

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)