

Product Name :
Manual Suction Unit (Foot / Pedal Operated)

Product Code :
NMSUC30HV3



Description :

Manual Suction Unit (Foot / Pedal Operated) Capacity: Foot Operated Vacuum pump giving vacuum more than 550mm (+ 50mm). Jars: Two Autoclavable PC jars with cap of 1 litre capacity. A suction cannula is also provided with the unit. The unit is fitted on a powder coated base

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,