

**Product Name :**  
Manual Suction Unit (Foot/Pedal Operated) - SUPERIOR

**Product Code :**  
NMSUC30HV4



**Description :**

Manual Suction Unit (Foot/Pedal Operated) - SUPERIOR Housing: Plastic Moulded Capacity: -600 mm Hg  $\pm$  10 , 200 ml / Stroke Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Dimension & Weight: 33 x 18 x 30 cms, 2.5 Kg. Overflow Safety: Mechanical Accessories: Tool Kit, Filter Paper & Push-in Nozzle

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,