

Product Name :
Suction Unit HI-VAC Plus (60 Ltrs./Min)

Product Code :
NMSUM10HV2



Description :

Suction Unit HI-VAC Plus (60 Ltrs./Min) Housing: Plastic moulded on two wheels Capacity: -600 mm Hg \pm 10 , 60 ltrs/ min Pump Type: Double Diaphragm Pump Jars: 2 x 2.5 Ltrs wide mouthed (Polycarbonate) with changeover lever. Provided with handle for hygiene while emptying/cleaning. Filter: Bacterial Filter Autoclavable/Reusable Tubing: 10 mm ID x 2 mtr. (PVC) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration Power: 220 V AC, 50 / 60 Hz, 180Watt (110V on request) Noise Level: 60 dB A \pm 3 Dimension & Weight: 46 x 34 x 86 cms, 22 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Foot Switch

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,