



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Hot Water Bottle Cover - Knitted, 100% Cotton

**Product Code :**  
MRG20100030



**Description :**

Hot Water Bottle Cover - Knitted, 100% Cotton Knitted Hot Water Bottle Cover, 100% Cotton Available in any color Washing Instructions: - Machine wash - Do not bleach - Tumble dry - Cool iron if required



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,