



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hot Water Bottle Cover - Knitted, Lambswool/Nylon

Product Code :
MRG20100031



Description :

Hot Water Bottle Cover - Knitted, Lambswool/Nylon Knitted Hot Water Bottle Cover, 80/20 Lambswool/ Nylon
Available in any color Washing Instructions: - Dry Clean Only



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,