



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Tourniquet - Latex

**Product Code :**  
MRG20100057



**Description :**

Tourniquet - Latex Features: Tourniquet Latex Elastic Strip Tubular Reusable, Non-sterile. Length: 75-100cm Approx. Instructions for use: A tourniquet is used to facilitate the puncture of veins. Safety Process: This item must be cleaned and disinfected as often as necessary



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,