



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :

Radiant Warmer Stand / Radiant Overhead Lamp
(Standard) with Microprocessor based temperature controller
with 2 modes Skin / Manual model

Product Code :

NMCBS26026



Description :

Radiant Warmer Stand / Radiant Overhead Lamp (Standard) with Microprocessor based temperature controller with 2 modes Skin / Manual model This is mobile Warmer Stand mounted on wheels. Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650Watts. Metallic overhead heater box. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards alogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures Essential Attachments Lower Shelf. S.S. IV Stand



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,