

**Product Name :**

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 3 modes Skin/Air/Manual

**Product Code :**

NMCBS26033

**Description :**

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 3 modes Skin/Air/Manual This is mobile Warmer Stand mounted on wheels with detachable infant care trolley

**Technical Specification**

**Heater Box** Heating Element : Ceramic Heater 650W with parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures

**Infant Care Trolley** Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder

**Mattress** density approx 21-25Kg/m<sup>3</sup> and with removable, washable, waterproof cover

**Acrylic X-ray tray** with X-ray guide for positioning the cassette

**2 Lower shelves** Mounted on wheel (100mm) The design aspects also allow easy dismantling baby care area during cleaning or fumigation

**Essential Attachments** Side Mounted Monitor tray, Pole for Syringe Pump / Infusion Pump S.S. IV Stand. 4 swivel castors (100mm) - 2 castors with brake