

Product Name :

Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode

Product Code :

NMCBS26052

**Description :**

Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode Technical Specification: Acrylic Single Wall Canopy with six port holes. Front Loading. Four Inlets for administration of O₂, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,