



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Phototherapy Stand

Product Code :
NMCBS26055



Description :

Phototherapy Stand Mobile Stand with height adjustment facility along with easy swiveling of source box. Unit can be used along with Infant Care Trolley, Radiant Warmer and Incubator 4 Special Blue CFL Tubes (20W each) 2 Daylight CFL Tubes (20W each) Single Head, surface Size, approx. 0.5m x 0.33m Head Height Adjustment, approx. 0.95 to 1.5m High irradiance at Skin level up to 40uw/cm²/nm (variable with distance) Wavelength - 420 to 500nm, with highest intensity at 470 nm Dual digital cumulative hour timer for tubes usage and patient exposure Antistatic Castors, 2 with brakes Power requirement – 230V/50Hz Power consumption – 100W



Medical Hospital Equipments Manufacturers

Website: SITE